ÖZET

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 LOT Bilgisi

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SUMMARY

MR-Detergent Solution is used to clean the electrodes flowcell as well as the pipeline of the ISE module on Mindray BS series automated biochemistry analysers.

Use the MR-Detergent Solution as instructed by the operators manual of the analyser.

APPEARANCE AND COMPONENTS

MR-Detergent Solution should be free of any precipitations, suspensions and floccules of any kind. It should be a slight yellow, transparent liquid. Do not use the product if it does not fit these specifications.

MR-Detergent Solution contains:
%3 ~ %5 w/w Sodium hypochlorite (effective chlorine concentration)

STORAGE AND STABILITY

MR-Detergent Solution is stable until its expiration date when its stored unopened at 5°C ~ 25 °C. The detergent vial should be kept away from light and should be capped tightly after use. Opened products are stable for 2 months if the product is kept under 5°C ~ 25 °C away from light.

LIMITATIONS

This product only applies to indirect ISE modules of Mindray BS series biochemistry analysers. In addition, the warranted results only come from processing tests strictly as instructed in the instrument manual and storing the product correctly.

WARNINGS

1. This product should not be used past the expiration date.
2. Wear proper personal protective equipment (such as gloves, safety goggles, lab coat etc.) and follow safe laboratory procedures when handling this product.
3. This product must be handled by skilled/trained medical professionals only.
4. Do not inhale or/and ingest this product. In case of inhaling or ingestion, consult a physician immediately. Avoid skin contact. In case of skin contact, flush the affected area with plenty of water. Avoid eye contact. In case of eye contact, flush the affected area with plenty of water and consult a physician immediately.
5. Dispose of any discarded material in accordance with the local government regulations.
6. Material safety data sheets are available upon request.

MR-Detergent Solution

APPLICATIONS

MR-Detergent Solution is only compatible with Mindrays BS series biochemistry analysers.

PRECAUTIONS

MR-Detergent Solution should only be used by trained medical professionals or laboratory staff for *in vitro diagnostic* use only.

HOW TO USE THIS PRODUCT

This product should be used as instructed on the user manual of your BS series biochemistry analyser.

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BİO AKTİF KİMYA İTH. İHR. İML. SAN. ve TİC. LTD. ŞTİ.